Sperm Washing Medium

with HSA

Catalog No. 9983

100 mL, 12 x 12 mL/Box

For sperm washing procedures.

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За postopke spiranja semenčic.

REFERENCES

Quinn, P. et. al.: Successful human in-vitro fertilization using a modified human tubal fluid medium lacking glucose and phosphate ions. Fert Steril 1995; 63: 922-24.

ENGLISH

EU Caution: For Professional Use Only

INDICATION FOR USE

Sperm Washing Medium with HSA is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of this medium for sperm washing.

DEVICE DESCRIPTION

This medium is a modification of Modified Human Tubal Fluid supplemented with Human Serum Albumin (HSA) (5 mg/mL) and containing the following components:

COMPOSITION:

Salts and Ions	Energy Sources
Sodium Chloride	Glucose
Potassium Chloride	Sodium Pyruvate
Magnesium Sulfate	Sodium Lactate
Potassium Phosphate	
Calcium Chloride	pH Indicator
	Phenol Red
Buffers	Protein Source
HEPES	Human Serum Albumin
Sodium Bicarbonate	
	Water
	WFI Quality

QUALITY ASSURANCE

Sperm Washing Medium is a culture medium which is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of Sperm Washing Medium is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

Sperm Washing Medium uses a buffering system composed of a 21 mM HEPES (N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid) and 4 mM Sodium Bicarbonate combination. This buffering system provides optimum pH maintenance over the physiologic range (7.2 to 7.4) and does not require the use of a CO₂ incubator.

PROTEIN SUPPLEMENTATION

Sperm Washing Medium contains 5 mg/mL human serum albumin which has been found to be non-reactive for Hepatitis B Surface Antigen (HB_sAg) and antibodies to Hepatitis C (HCV) and Human Immunodeficiency Virus (HIV) by approved testing methods. Donors of the source material have also been screened for Creutzfeldt-Jakob Disease (CJD)*.

DIRECTIONS FOR USE

The following is a general procedure for washing sperm from its surrounding seminal fluid.

- Bring medium to room temperature or 37°C.

(NOTE: Sperm Washing Medium should be tightly capped if warmed in a CO₂ incubator to avoid pH levels of 7.0 or less.)
- Allow the semen to liquefy at room temperature for 20 to 30 minutes.
- Using aseptic techniques, transfer the liquefied semen to a sterile 10 mL conical centrifuge tube and add 2 to 3 volumes of room temperature Sperm Washing Medium (for example, a 2 mL semen sample requires 4 to 6 mL of medium). Should the volume of the sperm medium mixture be greater than 5 mL, divide into two sterile tubes. By minimizing the volume per tube to 4 - 6 mL, the recovery of sperm will be maximized.

Samples having high viscosity may require a further processing to ensure total sperm recovery. (See Special Processing Considerations).

- Centrifuge the tubes at ambient temperature for 10 minutes using a g force of 200 - 300 x g.
- Using a sterile pipette, remove and discard the supernatant above the "sperm pellet" by aspiration. The sperm should then be resuspended by gently flicking the tube externally with the index finger. (Note: Do not use a vortex mixer for this step). Resuspend the sperm in 1 to 2 mL of fresh medium, recap and gently mix by inversion. Samples which were fractionated for the first centrifugation step should now be recombined into one tube.
- Recentrifuge as in Step 4.
- Using a sterile pipette, remove and discard the supernatant and resuspend the sperm pellet gently by manual agitation. Add fresh medium to a final volume of 0.5 mL. The sperm are ready for assisted reproductive procedures. (Note: The total volume of the nongravid uterus is 15 - 56 mL).

SPECIAL PROCESSING CONSIDERATIONS

Processing the highly viscous sample:

Some samples are naturally highly viscous even after liquefaction. These samples have the consistency of heavy syrup and may be among the most difficult to process.

- After the medium is added to an ejaculate, aspirate and expel the mixture gently using an 18 gauge needle and syringe. This will "shear" some of the viscous mucous.
- Limit the amount of medium-sperm mixture from Step 1 to 5 mL per centrifuge tube for the first centrifugation step.
- If after preprocessing the sample with the needle and syringe (Step 1), the sperm do not "pellet" in a normal manner (the sperm will appear as a "cloudy fiber" attached to the bottom of the centrifuge tube), carefully aspirate as much of the supernatant as possible without disrupting the "cloudy sperm fiber" using a sterile needle and syringe. This can be done by keeping the beveled edge of the needle firmly against the wall of the centrifuge tube and slowly start aspiration from the top of the tube downward. When as much of the supernatant as possible has been removed, add 2 or 3 mL of fresh medium. Repeat the process of drawing the mixture through the 18 gauge needle and syringe. Recentrifuge the mixture. The sperm should pellet normally after the second processing.
- On subsequent sample collections, the patient should be requested to produce a split ejaculate which will minimize the viscosity in the sperm rich portion of the sample.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles or vials refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

When stored as directed, Sperm Washing Medium is stable until the expiration date shown on the bottle label.

Duration Following Bottle Opening:
Product should be used within (8) weeks from opening when stored under the recommended conditions of 2° to 8°C.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle or vial of medium which shows evidence of particulate matter, cloudiness or is not pale rose in color.

To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that remains in the container after the sperm processing procedure is completed.

Do not use any bottle in which the sterile packaging has been compromised.

Sperm Washing Medium does not contain antibiotics. Antibiotic supplementation may be added just prior to use.

Appropriate precautions should be taken to ensure that the patient is not sensitized to the antibiotic of choice.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopeia specifications by established processes. It is strongly recommended that every time this product is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

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