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
Za postopke asistirane reprodukcije.

Glossary of Symbols:

	Catalog Number
	Lot Number
	Sterilized using aseptic processing techniques (filtration)
	Expiration: Year - Month - Day
	Caution, consult accompanying documents
	Consult instructions for use
	Storage Temperature 2-8°C
	Do Not Re-Sterilize
	Do Not Use If Package Is Damaged
	Phthalate, DBP, DEHP
	Manufacturer
	U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CE Mark
	Emergo Europe - Westervoortsedijk 60 6827 AT Arnhem The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

 FUJIFILM Irvine Scientific, Inc.

2511 Daimler Street, Santa Ana, California 92705 USA

Telephone: 1 949 261 7800 • 1 800 437 5706 • Fax: 1 949 261 6522 • www.irvinesci.com

© 2023 FUJIFILM Irvine Scientific, Inc. All rights reserved. FUJIFILM Irvine Scientific, and its logo, Vit Kit, and CryoTip, are registered trademarks of FUJIFILM Irvine Scientific, Inc. in various jurisdictions. All other trademarks are the property of their respective owners.

PN 40696 Rev.23

Effective Date: 31-JUL-2023

CryoTip Diagrams:

Figure 1:


Figure 2:


Figure 3:




Figure 4:

TS = Thawing Solution
DS = Dilution Solution
WS = Washing Solution

HSV Device Diagrams:

Cryolock Diagrams:

Figure 1:



Figure 2:

ENGLISH

EU CAUTION: For Professional Use Only.

INDICATION FOR USE

Vit Kit-Thaw (Vitrification Thaw Kit) is intended for use in the thawing of vitrified oocytes (MII), pronuclear (PN) zygotes through day 3 cleavage stage embryos and blastocyst stage embryos that have been vitrified using Vitrification Freeze Kit (Catalog #90133-SO)

DEVICE DESCRIPTION

Thawing Solution-TS is a HEPES buffered solution of Medium-199 containing Gentamicin Sulfate, 1.0 M sucrose and 20% (v/v) Dextran Serum Supplement (DSS).

Dilution Solution-DS is a HEPES buffered solution of Medium-199 containing Gentamicin Sulfate, 0.5M sucrose and 20% (v/v) DSS.

Washing Solution-WS is a HEPES buffered solution of Medium-199 containing Gentamicin Sulfate and 20% (v/v) DSS.

DSS is a protein supplement consisting of 50 mg/mL therapeutic grade Human Serum Albumin (HSA) and 20 mg/mL Dextran. DSS is used at 20% (v/v) in Vit Kit-Thaw for a final concentration of 10 mg/mL HSA and 4 mg/mL Dextran.

These three solutions are to be used in sequence according to the step-wise microdrop warming protocol.

COMPOSITION

<u>Salts and Ions</u>	Nicotinic Acid Amide	<u>Protein Source</u>	Lysine	Hydroxyproline
Sodium Chloride	Pantothenic Acid	Human Serum Albumin	Proline	Cystine
Sodium Phosphate	Riboflavin	<u>Buffers</u>	Tyrosine	<u>Other</u>
Potassium Chloride	Pyridoxine	Sodium Bicarbonate	Alanine	Guanine
Magnesium Sulfate	Thiamine	HEPES	Aspartic Acid	Hypoxanthine
Sodium Acetate	Biotin	<u>pH Indicator</u>	Glutamic Acid	Thymine
Calcium Chloride	Alpha-Tocopherol	Phenol Red	Isoleucine	Uracil
Ferric Nitrate	Sodium Bisulfite	<u>Macromolecules</u>	Leucine	Xanthine
Choline Chloride	<u>Antioxidant</u>	Sucrose	Methionine	Adenosine
<u>Vitamins and Minerals</u>	Glutathione	Dextran	Phenylalanine	Adenine Sulfate
Ascorbic Acid	<u>Antibiotic</u>	<u>Amino Acids</u>	Serine	Deoxyribose
Aminobenzoic Acid	Gentamicin Sulfate	Arginine	Threonine	Ribose
Calciferol	<u>Energy Substrates</u>	Glycine	Tryptophan	<u>Water</u>
Folic Acid	Glucose	Histidine	Valine	WFI Quality
Nicotinic Acid	Inositol		Cysteine	

QUALITY ASSURANCE

The solutions in Vit Kit-Thaw are membrane filtered and aseptically processed according to manufacturing procedures which have been validated.

Each lot of Vit Kit-Thaw receives the following tests:

Solutions:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology (≤ 0.6 EU/ml)
- Mouse Embryo Assay (one-cell) ($\geq 80\%$ expanded Blastocyst)
- Sterility by the current USP Sterility Test <71> (Pass)

All results are reported on a lot specific Certificate of Analysis which is available upon request.

FOR WARMING CRYOTIP AS THE CARRIER:

MATERIALS REQUIRED BUT NOT INCLUDED

- Connector (Catalog #40736) or adaptor
- Sterile petri dishes (50 X 9 mm, Falcon 351006 or equivalent)
- Disposable gloves
- Hamilton GASTIGHT® Syringe (50 μ l) catalog #80901
- Transfer pipettes (pulled glass pipettes or micropipette tips with an inner tip diameter of ~ 200 μ m)
- Tweezers
- Stopwatch or timer
- Liquid nitrogen reservoir (dewar or styrofoam container with lid, 1-2 L volume)
- Liquid nitrogen (sufficient volume to achieve 4 inch depth in reservoir)
- Sharp scissors (sterile)
- 37°C waterbath
- Culture medium with protein, pre-equilibrated to 37°C in CO₂ incubator prior to thawing procedure.
- 37°C incubator without CO₂, or heating stage.

DIRECTIONS FOR USE

Vit Kit-Thaw components (per application):

- 50 μ l of Thawing Solution –TS
- 50 μ l of Dilution Solution-DS
- 100 μ l of Washing Solution-WS
- 1 Connector

WARMING PROTOCOL

(FOR OOCYTES AND EMBRYOS):

NOTE: Procedures are to be done at room temperature (20-27°C). DO NOT use heated microscope stage for the following procedures.
CAUTION: Minimize exposure of specimen to light during manipulations through thawing solutions.

1. Bring the quantity to be used of TS, DS and WS to room temperature (20-27°C) prior to warming vitrified specimens.
NOTE: Avoid bringing the entire vials of TS, DS and WS to room temperature repeatedly when a small quantity of the solution is needed each time. It is better to aliquot the quantity to be used and return the vials to 2-8°C right after aliquoting.
2. Fill the liquid nitrogen reservoir with liquid nitrogen (~80 % full) and place close to the LN₂ freezer containing the specimens to be thawed.
3. Remove the canes with goblets containing the CryoTips with vitrified specimens from liquid nitrogen storage and transfer them into the liquid nitrogen filled reservoir.
CAUTION: Make sure that CryoTips remain submerged in LN₂ (in goblet) during transfer from storage to LN₂ reservoir to prevent uncontrolled thawing of specimens.
Place the reservoir close to the microscope for rapid manipulation.
4. Label a sterile petri dish (or lid) with necessary information.
5. Gently invert each vial of TS, DS and WS twice to mix contents before use.
6. Prepare the dish with droplets of solutions for warming protocol as follows:
Aseptically dispense a sequence of 2 microdrops on an inverted lid of sterile Petri dish as shown in Figure 1, and place the dish on the microscope stage:
 - one-50 µl drop of TS
 - one-50 µl drop of DS
 - (Two drops of WS will be set up later at step 11)
7. Place the 37°C waterbath close to the microscope. Have the following nearby: a transfer pipet and tips, sterile sharp scissors, Hamilton syringe and sterile wipes.
8. Using tweezers (or tongs), retrieve the specific CryoTip from the cane in liquid nitrogen, quickly immerse the CryoTip into the 37°C waterbath (≥ 500 ml) and swirl gently for 3 seconds to warm (see Figure 2) at + 24,000°C/min.
9. Quickly dispense contents of CryoTip as follows (see Figure 3):
 - Quickly wipe the CryoTip dry with a sterile wipe
 - Remove metal cover sleeve
 - Cut the seal on the wide end of the CryoTip at Mark #4
 - Attach wide end of the CryoTip securely to Hamilton syringe or appropriate aspiration tool using a connector or adaptor.
NOTE: Lift up the plunger of the syringe about 0.5 inches before attaching the syringe to the Connector and the CryoTip.
 - Gently wipe the fine tip dry with a sterile wipe.
 - With the CryoTip positioned over the prepared thawing dish, quickly cut the seal on the Mark #2 of the fine tip end and dispense the contents of the CryoTip as a small drop (~ 1 µl) onto a dry region of the dish above the TS drop (see Figure 4).
NOTE: AVOID BUBBLES WHILE DISPENSING THE CONTENTS.
10. Merge the TS drop with the CryoTip content and allow gradual mixing for 1 minute (see Figure 4).
Note: the specimens will shrink and float to the top of the drop.
NOTE: After each transfer of the specimen(s), blow out any remaining fluid in the transfer pipet and draw up some solution from the next drop prior to the next manipulation. Avoid creating bubbles during the transfers.
11. Draw up some DS into the transfer pipet and transfer the specimen(s) from the drop of TS with minimal volume to the drop of DS for 4 minutes.
NOTE: The specimen will remain shrunken during exposure to DS.)
DURING THIS TIME SET UP THE TWO-50 µl DROPS OF WS (WS1, WS2), AS SHOWN IN FIGURE 4.
12. Transfer the specimen(s) to the drop of WS (WS1) for 4 minutes.
NOTE: The specimen(s) should re-expand to the original size within 2-3 minutes in WS.
13. Then transfer the specimen(s) to the second drop of WS (WS2) for 4 minutes.
14. Transfer warmed OOCYTE(S) to pre-equilibrated culture medium with 20% (v/v) protein supplement or 12 mg/ml for recovery (2-3 hours to allow time for spindle re-formation) prior to subsequent manipulations.
There are two options for warmed EMBRYO(S):
 - a) For immediate transfer to patient: transfer EMBRYO(S) to pre-equilibrated 'transfer' medium containing 20% (v/v) protein supplement or 12 mg/ml.
 - b) For further culture: transfer EMBRYO(S) to pre-equilibrated culture medium containing 20% (v/v) protein supplement or 12 mg/ml for a 4 hour recovery period. After recovery period, transfer EMBRYO(S) to culture medium with 10% (v/v) protein and incubate accordingly until desired developmental stage has been reached for transfer to patient.

FOR WARMING HSV DEVICE AS THE CARRIER:

MATERIAL REQUIRED BUT NOT INCLUDED

- Sterile 4-well dish (Nunc 179830, 144444 or equivalent), or organ culture dish (BD Falcon 353037)
- Disposable gloves
- Transfer pipettes
- Tweezers
- Stopwatch or timer
- Liquid nitrogen reservoir
- Liquid nitrogen
- Scissors, Knipex or other wire cutter device
- Culture medium with protein, pre-equilibrated to 37°C in CO₂ incubator prior to thawing procedure.
- 37°C incubator without CO₂, or heating stage

DIRECTIONS FOR USE

Vit Kit-Thaw components (per application)

- 250 µl of Thawing Solution –TS
- 50 µl of Dilution Solution – DS
- 100 µl of Washing Solution – WS

WARMING PROTOCOL

(FOR OOCYTES AND EMBRYOS):

NOTE: The warming steps include plunging the device into the 37°C TS and subsequent diluting and washing in DS and WS at room temperature

1. Set up warming dishes (as shown in HSV device diagram):
 - At 37°C: Aseptically dispense 250 µl of TS into a sterile 4-well dish or an organ culture dish and place it in a 37°C incubator without CO₂ or on a heating stage at least 30 minutes prior to warming procedure

NOTE: For oocytes, dispense a minimum of 1 ml of TS

2. Identify HSV Straw(s) to be warmed from LN₂ storage and quickly transfer to LN₂ filled reservoir in preparation for thawing procedure.
3. Place LN₂ reservoir close to microscope for subsequent rapid manipulation.
4. Remove TS dish from 37°C incubator or heating stage and place it under focus on top of the microscope stage.
5. Lift the straw enough to expose the colored handling rod. Make sure the end with the specimen(s) remains immersed in the LN₂.
6. Use a Knipex (or other wire cutter device) to cut the straw at the height of the colored handling rod. The red cut-length guide on the Knipex should be positioned in maximum length position or removed.
 - Alternatively, use fingers and thumb to spin the straw while making cutting movements with scissors, 10 mm under the top of the colored handling rod.
7. With one swift but controlled motion, quickly grab the handling rod and extract it out of the straw.
8. Immediately plunge the gutter into the 37°C TS and gently swirl to detach specimens from device and leave it for 1 minute.

Steps 9-12 must be performed at room temperature (22-27°C).

- At room temperature: Aseptically dispense one (1) 50 µl drop of DS on a sterile Petri dish

9. Draw up some DS into the transfer pipette and transfer the specimen(s) from the TS drop with minimal volume to the drop of DS for 4 minutes.

NOTE: The specimen will remain shrunken during exposure to DS.

DURING THIS TIME SET UP THE TWO-50 µl DROPS OF WS (WS1, WS2), AS SHOWN IN HSV DEVICE DIAGRAM.

10. Transfer the specimen(s) to the drop of WS (WS1) for 4 minutes.

NOTE: The specimen(s) should re-expand to the original size within 2-3 minutes in WS.

11. Then transfer the specimen(s) to the second drop of WS (WS2) for 4 minutes.
12. Transfer warmed OOCYTE(S) to pre-equilibrated culture medium with 20% (v/v) protein supplement or 12 mg/ml for recovery (2-3 hours to allow time for spindle re-formation) prior to subsequent manipulations.

There are two options for warmed EMBRYO(S):

- a) For immediate transfer to patient: transfer EMBRYO(S) to pre-equilibrated 'transfer' medium containing 20% (v/v) protein supplement or 12 mg/ml.
- b) For further culture: transfer EMBRYO(S) to pre-equilibrated culture medium containing 20% (v/v) protein supplement or 12 mg/ml for a 4 hour recovery period. After recovery period, transfer EMBRYO(S) to culture medium with 10% (v/v) protein and incubate accordingly until desired developmental stage has been reached for transfer to patient.

NOTE: Recovered oocytes must be fertilized using ICSI for optimal fertilization after vitrification.

FOR WARMING CRYOLOCK™ AS THE CARRIER:

MATERIALS REQUIRED BUT NOT INCLUDED

- Sterile 4-well dish or sterile small petri dishes (35 X 10 mm or equivalent)
- Disposable gloves
- Transfer pipettes
- Stopwatch or timer
- Liquid nitrogen reservoir
- Liquid nitrogen
- Culture medium with protein, pre-equilibrated to 37°C in CO₂ incubator prior to thawing procedure
- 37°C incubator without CO₂ or heating stage
- Forceps

DIRECTIONS FOR USE

Vit Kit-Thaw components (per application)

- 1 ml of Thawing Solution – TS
- 50 µL of Dilution Solution – DS
- 100 µL of Washing Solution – WS

WARMING PROTOCOL

NOTE: The warming steps include plunging the device into the 37°C TS and subsequent diluting and washing in DS and WS at room temperature

1. Set up thawing dish (as shown in Cryolock diagram, Figure 1):
 - At 37°C: Aseptically dispense a minimum volume of 1 ml of TS and warm to 37°C in an incubator without CO₂ or on a heating stage at least 30 minutes prior to starting warming procedure.
2. Identify the Cryolock sample(s) to be warmed and quickly transfer from LN₂ storage to an LN₂ filled holding reservoir in preparation for warming procedure.
3. Place LN₂ filled holding reservoir in close proximity to the working area and stage of the microscope in order to achieve subsequent rapid manipulation from reservoir to TS.
4. Remove TS dish from 37°C incubator or heating stage and place it under focus on top of the microscope stage.
5. Using forceps, hold the upper end of the Cryolock body with the identification label facing up.
 - Option A: Quickly but gently remove the cap under LN₂, twisting the parts until release.
 - Option B: Quickly take the Cryolock out from the LN₂, then quickly remove the cap with a gentle twist.

NOTE: Laboratory should consult their own procedures and protocols. Option A is not cleared for use in the U.S.

6. Immediately plunge the concave tip of the Cryolock, with the specimen(s) facing up, into the 37°C TS. Under microscopic observation, gently move the Cryolock until the specimen(s) are released from the tip.
7. Leave the specimen(s) for a total of 1 minute in the TS.
8. Thirty (30) seconds following the initial plunge, gently pipette the specimen(s), if floating, and place at the bottom of the TS.

Steps 9-12 must be performed at room temperature (22-27°C).

- At room temperature: Aseptically dispense one (1) 50 µL drop of DS on a sterile petri dish (see Cryolock diagram, Figure 2)

9. Transfer specimen(s) to DS for 4 minutes. Gently pipette specimens once to ensure complete rinse in DS.

NOTE: The specimen will remain shrunken during exposure to DS.

10. During the 4 minute exposure in DS, aseptically dispense two (2) 50 µL drops of WS (WS1, WS2) as shown in diagram.

11. Transfer specimen(s) to WS1 then WS2 for 4 minutes each, undisturbed.

NOTE: The specimen(s) should re-expand to the original size within 2-3 minutes in WS.

12. Transfer warmed OOCYTE(S) to pre-equilibrated culture medium with 20% (v/v) protein supplement or 12 mg/ml for recovery (2-3 hours to allow time for spindle re-formation) prior to subsequent manipulations.

There are two options for warmed EMBRYO(S):

- a) For immediate transfer to patient: transfer EMBRYO(S) to pre-equilibrated 'transfer' medium containing 20% (v/v) protein supplement or 12 mg/ml.
- b) For further culture: transfer EMBRYO(S) to pre-equilibrated culture medium containing 20% (v/v) protein supplement or 12 mg/ml for a 4 hour recovery period. After recovery period, transfer EMBRYO(S) to culture medium with 10% (v/v) protein and incubate accordingly until desired developmental stage has been reached for transfer to patient.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened vials refrigerated at 2°C to 8°C. When stored as directed, Vitrification Thaw Kit Solutions are stable until the expiration date shown on the vial labels.

Do not use media for more than eight (8) weeks once containers have been opened.

As human source material is present in the product it may develop some particulate matter during storage. This type of particulate matter is not known to have an effect on product performance.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any vial of solution which shows evidence of damage, leaking, particulate matter, cloudiness or has shifted color. Discard the product in accordance with applicable regulations.

To avoid problems with contamination, handle using aseptic techniques.

Currently, research literature indicates the long-term effects of vitrification on oocytes and embryos remains unknown.

Do not use any bottle in which the sterile packaging has been compromised.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: This product contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

CONTRAINDICATION

Product contains Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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