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Za postupke asistirane reprodukcije.

REFERENCES

Biggers, JD, and Racowsky, C. The development of fertilized human ova to the blastocyst stage in KSOM^{AA} medium: is a two-step protocol necessary? RBMOnline, 5:133-140, 2002.

Pool, TB. Recent advances in the production of viable human embryos in vitro. RBMOnline, 4:294-302, 2002.

Biggers, JD. Thoughts on embryo culture conditions. RBMOnline, 4 (suppl.1):30-38, 2001.

Lane, M., Hooper, K., and Gardner, DK. Effect of essential amino acids on mouse embryo viability and ammonium production. J. Asst. Reprod. Genet. 18: 519-525, 2001

Biggers, JD, and McGinnis, LK. Evidence that glucose is not always an inhibitor of mouse preimplantation development in vitro. Hum. Reprod 16:153-163, 2001.

Devreker, F., Van den Bergh, M., Biramane, J., Winston, RML., Englert, Y., and Hardy, K. Effects of taurine on human embryo development in vitro. Hum. Reprod. 14: 2350-2356, 1999.

REF

LOT

STERILE A

Catalog Number

Lot Number

Sterilized using aseptic processing techniques (filtration)

Expiration:
Year - Month - Day

Caution, consult accompanying documents

Consult instructions for use

Storage Temperature
2-8°C

Manufacturer

Do not resterilize

Do not use if package is damaged

U.S. Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

CE Mark

Emergo Europe - Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

*Symbol Reference - EN ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labeling.

ENGLISH

EU CAUTION: For Professional Use Only

INDICATION FOR USE

Continuous Single Culture Complete (CSCM-C) is intended for use in assisted reproductive procedures which include human gamete and embryo manipulation. These procedures include the use of CSCM-C as a culture medium from fertilization through day 5/6 of development in vitro.

DEVICE DESCRIPTION

Continuous Single Culture Complete is a novel solution optimized for use in an uninterrupted culture system without dish change or medium renewal.

CSCM-C contains Human Serum Albumin (HSA) for a final total protein concentration of 5 mg/mL and the antibiotic Gentamicin Sulfate (10 µg/mL).

COMPOSITION:	
Salts & Ions	Amino Acids
Sodium Chloride	Alanine
Potassium Chloride	Asparagine
Calcium Chloride	Aspartic Acid
Magnesium Sulfate	Glutamic Acid
Potassium Phosphate	Glycine
Buffer	Proline
Sodium Bicarbonate	Serine
	Arginine
Energy Substrates	Cystine
Sodium Pyruvate	Histidine
Glucose	Isoleucine
Sodium Lactate	Leucine
Antioxidant	Lysine
EDTA	Methionine
Sodium Citrate	Phenylalanine
Dipeptide	Threonine
Alanyl-Glutamine	Tryptophan
pH Indicator	Tyrosine
Phenol Red	Protein
Antibiotic	Human Serum Albumin
Gentamicin Sulfate	Water
	WFI Quality

QUALITY ASSURANCE

CSCM-C is membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of CSCM-C is tested for:

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology (≤ 0.25 EU/mL)
- Biocompatibility by Mouse Embryo Assay (one-cell at ≥80% expanded blast 96h)
- Sterility by the current USP Sterility Test <71> Human Sperm Survival Assay (HSSA) (≥70% motility at 24h).

All results are reported on a lot specific Certificate of Analysis which is available upon request.

BUFFER SYSTEM

CSCM-C uses sodium bicarbonate as a buffering system. This is specifically designed for use in a CO₂ incubator.

DIRECTIONS FOR USE

CSCM-C is a complete, ready-to-use medium containing Human Serum Albumin (HSA). It is not necessary to add any protein before use.

EQUILIBRATION

CSCM-C should be pre-warmed to 37°C and equilibrated to the desired pH overnight in a 5-6% CO₂ incubator prior to use. A sufficient volume of medium is required so that oocyte recovery, insemination and embryo culture dishes can be prepared.

The following are general procedures for the indications for use of CSCM-C.

Fertilization:

On the day before oocyte retrieval, prepare insemination dishes with CSCM-C overlaid with oil and pre-equilibrate overnight to 37°C in a CO₂ incubator. Immediately upon oocyte collection and identification, place oocytes into the pre-equilibrated medium and return to the incubator for the desired period (1-4 hours) prior to insemination by conventional IVF or ICSI.

- Conventional IVF:**
- It is recommended to aseptically dispense 50,000 - 100,000/mL motile sperm per microdroplet containing 1-3 oocytes.
 - Return the insemination dish to the incubator and check for normal fertilization 16-20 hours post insemination.
- Intracytoplasmic Sperm Injection (ICSI):**
- Following at least 1 hour post oocyte denuding (and no more than 4 hours following oocyte retrieval), remove denuded oocytes from incubator and inseminate with sperm per standard ICSI protocol for your individual laboratory.
 - Immediately following insemination, place 1-3 inseminated oocytes into a fresh drop of the pre-equilibrated insemination dish, return dish to the incubator and check for normal fertilization 16-20 hours post insemination.

Embryo Culture:

On the day of fertilization (one day prior to fertilization assessment), prepare embryo culture dishes with CSCM-C overlaid with oil and pre-equilibrate overnight to 37°C in a 5-6% CO₂ incubator.

Following fertilization assessments with the identification of the presence of normal fertilization (two pronuclei and two polar bodies), transfer 2PN zygotes into the pre-equilibrated CSCM-C culture dish previously prepared. It is recommended to allow the embryos to grow in a continuous, uninterrupted culture system, without changing medium, until the desired developmental stage is reached (up to day 5/6 of development).

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened bottles refrigerated at 2° to 8°C.

Do not freeze or expose to temperatures greater than 39°C.

Duration Following Bottle Opening:

Product should be used within eight (8) weeks from opening.

As human source material is present in the product it may develop some protein precipitate during storage or in culture over time. This protein precipitate is not known to have an effect on product performance.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures. These procedures include the intended application for which this device is intended.

The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Do not use any bottle of medium which shows evidence of particulate matter, cloudiness or is not pale orange in color.

Not for injection use.

CAUTION: To avoid problems with contamination, handle using aseptic techniques and discard any excess medium that shows any evidence of contamination after opening.

EU: Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. It is strongly recommended that every time FUJIFILM Irvine Scientific, Inc. Reproductive Media Products culture media are administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

US: CSCM-C contains Human Serum Albumin (HSA). Human source material used in the manufacture of this product has been tested by FDA-licensed kits and found to be non-reactive to the antibodies to Hepatitis C (HCV), and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all human source material as if it were capable of transmitting infection, using universal pre-cautions. Donors of the source material have also been screened for CJD.

CONTRAINDICATION

CSCM-C contains the antibiotic Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

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