


MHM™ :


A Unique and Improved IVF Handling Media

In Vitro Stressors



Reduce stress to improve embryo development and ART outcomes

What is pH?




- Acids increase the concentration of hydrogen ions
- Bases decrease the concentration of hydrogen ions

pH is the measure of $[H^+]$

Internal pH (pHi)

- Cells contain pHi regulatory mechanisms
 - $\text{HCO}_3^-/\text{Cl}^-$ exchanger $>7.2-7.3$
 - Na^+/H^+ antiporter <6.8
 - Na^+ dependent $\text{HCO}_3^-/\text{Cl}^-$ exchanger <7.0
- pHi follows external pH (pHe) of media initially**




pHi and the Embryo

- Slightly raising or lowering pHi for 3hrs results in disorganization of mitochondria and actin cytoskeletal elements (Squirrell et al. 2001)
 - Regulate development and chromosome dynamics
- Raising pHi 0.09-0.15 for 4hrs significantly changes metabolism (Lane et al. 2000)
 - Metabolism is correlated with developmental competence
- Lowering pHi ~0.15 affects blastocyst development and resulting fetal size (Zander-Fox et al. 2010)

pH and ART

- Denuded mature oocytes lack robust pH_i regulatory mechanisms
 - Activated ~6h after fertilization

(Phillips et al. 1998, 2000, 2002)




- Cryopreserved/thawed embryos have reduced ability to regulate pH_i

- ~3h recovery (Lane et al. 2000)




- Sperm pH_i and function are influenced by pH_e (Hamamah et al 1996)



Proper and stable pH_e is crucial

External Media pH (pHe)


Incubator vs. Media



How do we control pHe outside the incubator?

IVF Handling Media

Media that uses reduced bicarbonate concentration and includes a zwitterionic buffer, like HEPES or MOPS, to maintain pH outside the incubator



Importance of Handling Media

- Brief exposure to inappropriate handling media can significantly reduced embryo development
 - Hamster (Escriba et la. 2001)
 - Rabbit (Farrell & Bavister 1984)
 - Cow (Palasz et al. 2008)
 - Mouse (Gardner & Lane 1996)
 - Human (Morgia et al. 2006)

Common Concerns with Buffers

- Buffers, like HEPES, are toxic
- Injection of buffers may alter pH (Morgia et al. 2006)
- HEPES and MOPS block Cl⁻ channels and may inhibit blastocyst development (Yamamoto and Suzuki, 1987, Butler et al., 1988).
- Cell specific sensitivity to particular buffers (Eagle 1971)
- Concentration dependent side-effects of buffers (Downs & Mastropoki 1997, Iwasaki et al. 1999)

Though many of the concerns are unfounded, this presents an opportunity to develop an improved handling medium

Objective

- Develop a unique and improved IVF handling medium
 - Accomplish with minor modifications to an already accepted medium to facilitate acceptance
 - 1) Reduced buffer concentration
 - 2) Improved buffer selection
 - Buffering capacity (pKa)
 - 3) Inclusion of select amino acids

Reduced Buffer Concentration

Concerns with Buffers

- Detrimental effects described with some buffers may be concentration dependent (Downs & Mastropoki 1997, Iwasaki et al. 1999)
 - Increasing buffer from 20 to 25mM prevented pharmacologic inhibition of oocyte maturation
 - Increasing buffer above 35mM increased pig embryo degeneration

Buffer Concentration

- ~2x the concentration of zwitterionic buffer as bicarbonate is sufficient to stabilize pHe
(Freshney 1983)
- Most IVF handling media contain ~21mM buffer
 - 5.25X conc. of 4mM bicarb


How did we arrive at current formulations?

Objective


Determine if reduced buffer concentration in IVF handling media maintains pH stability and supports embryo development

MHM™ - pH Stability

Dr. Cassuto's Lab




Dr. Swain's Lab



MHM™ - MEA Development

2 Beta Site Labs – 2 hr MHM Exposure – 96h of culture in incubator



Custom Combination Buffering System to Optimize Buffering Capacity Over a Range of Temperatures

Objective

Formulate a dual buffered system that offers improved buffering capacity over a range of temperatures compared to current single buffered media containing only HEPES or MOPS

Buffer Selection

- Buffers are selected based on ability to support cell growth
 - Not all buffers are compatible with all cell types
- Compatible buffers are then chosen based on their maximal buffering capacity...or ability to maintain a specific and stable pH
 - Maximal buffering is indicated by a buffer's pKa value
 - Maximal buffering is obtained when pKa is equal to the the desired pH (7.2-7.4 in IVF labs)


Maximal Buffering: $\text{pH} = \text{pK}_a$

Buffers & pK_a

| Common Name | pK _a at 20°C |
|-------------|----------------------------|
| TAPSO | 7.7 |
| DIPSO | 7.6 |
| HEPES | 7.55 |
| TES | 7.5 |
| Phosphate* | 7.21 |
| MOPS | 7.20 |
| Carbonate* | 6.38 |

Temperature in the lab

Surface Temperature 37.0°C



Temp will vary throughout the lab, and even depending on the dish/volume used

Buffers & pKa

Temperature Impacts Buffering

| Common Name | pK _a at 20°C | pK _a at 37°C |
|-------------|----------------------------|----------------------------|
| TAPSO | 7.7 | 7.39 |
| DIPSO | 7.6 | 7.35 |
| HEPES | 7.55 | 7.31 |
| TES | 7.5 | 7.16 |
| Phosphate* | 7.21 | 7.19 |
| MOPS | 7.20 | 6.95 |
| Carbonate* | 6.38 | 6.30 |

Article

New pH-buffering system for media utilized during gamete and embryo manipulations for assisted reproduction

Jason E Swain¹, Thomas B Pool

Fertility Center of San Antonio, San Antonio, TX 78229 USA


¹Correspondence: e-mail: swainj@med.umich.edu




Combine zwitterionic buffers

- To avoid concern with elevated concentration and possible toxicity:
- To optimize pH buffering capacity (pKa) considering temp

Combination Buffers



Combination pH Buffering



Conclusion

- Varying ratios of HEPES and MOPS allows for adjustment of optimal buffering and allows for the creation of a unique handling media that provides optimal buffering over the range of temperatures encountered in and IVF lab

Amino Acids and IVF Handling Media

Objective

- **Examine the effect of various amino acids and their combinations on embryo development when included in an IVF handling medium, focusing on amino acids that could be included in a universal handling media suitable for gametes and all embryo stages**

- 1) NEAA
- 2) Glutamine
- 3) Glycine
- 4) Taurine

Amino Acids

- Amino acids act as metabolic substrates, osmolytes and regulators of pH_i (Lane 2000)
- Absence of amino acids in handling media resulted in significantly decreased blastocyst formation in mouse (Gardner & Lane 1996)

All media should contain some assortment of amino acids

Amino Acids

- Some amino acids are beneficial, while others are detrimental – dependent upon concentration and cell type/stage
 - Cleavage stage mouse embryos benefit from inclusion of NEAA (glutamine), while EAA are beneficial post-compaction (Gardner & Lane 1993, 1994, 1997a,b)
 - Glycine, taurine and glutamine found most beneficial for hamster embryos (McKiernann et al. 1995)
 - Taurine acts as an osmolyte and is beneficial for human embryos (Dumoulin et al. 1997, Dawson & Baltz 1997)
 - 1 report of benefit of glutamine for human embryos grown in glucose-free media (Devreker et al. 1998)
 - Glycine is a potent osmolyte and transporter identified in human embryos (Hammer et al. 2000)

Amino Acids - Results

- After numerous experiments:
 - No significant benefit of NEAA alone or in combination with other amino acids were found at varying concentrations
 - No significant benefit was found when including glutamine, alone, or in combination with other amino acids
 - Though not significant, taurine and glycine supplemented in combination gave slightly higher rates of embryo development compared to other treatments
 - Are included in most IVF medium, have known functions

Why Exclude Glutamine?


- Glutamine impacts glucose metabolism (Chatot et al. 1990, Du & Wales 1993)
 - Trying to avoid metabolic perturbations – saw no benefit in our study
- Glutamine is labile in culture and can form harmful ammonia – this necessitates use of dipeptide forms
 - Dipeptides may not function as optimally as individual amino acids (Swain et al. 2011)
- Glutamine utilizes the same transporter as glycine for osmoregulation - redundant
 - Glycine has been shown to inhibit glutamine transport in post compaction mouse embryos, likely because both use the same GLYT1 transporter (Richards *et al.*, 2010)

Rationale for Amino Acid Selection

- NEAAs don't appear to be beneficial in the context of the basal media used
 - Why include unnecessary amino acids and risk potential negative side effects like ammonia buildup?
- No significant benefit of glutamine observed in context of our basal medium and potential drawbacks exist
- glycine + taurine appear slightly beneficial
 - Known/proven osmolytes and/or benefit in human embryos


MHM™ - Osmo Protection

■ mHTF 320mOsm ■ MHM™ 320mOsm




MHM™ Clinical Testing

MHM™ 1-cell MEA



MHM™ - Clinical ISCI Data

Human Sperm Motility following 24h Culture



MHM™ - Clinical ISCI Data

- Rationale for Testing with ICSI*** - Oocyte is most pHe sensitive cell stage
- Most invasive use of buffered media and most likely scenario to see impact

**3 IVF Clinics
47 Patients – 594 Oocytes**

